

## Rudolf Mask Replacements Cross Reference

Mask Number: Complete mask, headcap & valve	Description of Mask	Face Piece ONLY	Adapter: Facepiece to valve	Valve	Ports OD
<b>Adult Mouth-only Breathing Masks</b>					
R7921	Lg with lg Y-valve	R669040	R200542	R2730	35mm
R7922	Lg with lg T-valve	R669040	R200542	R2700	35mm
R7923	Lg with med Y-valve	R669040	R200543	R2630	35mm
R7924	Lg with med T-valve	R669040	R200543	R2600	35mm
R7931	Med with lg Y-valve	R669041	R200542	R2730	35mm
R7932	Med with lg T-valve	R669041	R200542	R2700	35mm
R7934	Med with med T-valve	R669041	R200543	R2600	35mm
R7936	Med with sm T-valve	R669041	R200591	R1410A	22mm
R7941	Sm with lg Y-valve	R669042	R200542	R2730	35mm
R7942	Sm with lg T-valve	R669042	R200542	R2700	35mm
R7944	Sm with med T-valve	R669042	R200543	R2600	35mm
R7946	Sm with sm T-valve	R669042	R200591	R1410A	22mm
V86018	Adult Headcap, straps & clips				
<b>Pediatric Mouth-only Breathing Masks</b>					
R7951	Lg ped with med Y-valve	R669044	R200584	R2630	35mm
R7952	Lg ped with med T-valve	R669044	R200584	R2600	35mm
R7953	Lg ped with sm Y-valve	R669044	R200585	R1420A	22mm
R7954	Lg ped with sm T-valve	R669044	R200585	R1410A	22mm
R7973	Sm ped with sm Y-valve	R669071	R200585	R1420A	22mm
R7974	Sm ped with sm T-valve	R669071	R200585	R1410A	22mm
R201079	Large ped Headcap, straps & clips				
R201133	Small ped Headcap, straps & clips				
<b>Adult Mouth and Nose Breathing Masks</b>					
R8921	Lg with lg Y-valve	R669050	R200542	R2730	35mm
R8922	Lg with lg T-valve	R669050	R200542	R2700	35mm
R8923	Lg with med Y-valve	R669050	R200543	R2630	35mm
R8924	Lg with med T-valve	R669050	R200543	R2600	35mm
R8931	Med with lg Y-valve	R669051	R200542	R2730	35mm
R8932	Med with lg T-valve	R669051	R200542	R2700	35mm
R8934	Med with med T-valve	R669051	R200543	R2600	35mm
R8936	Med with sm T-valve	R669051	R200543	R1410A	22mm
R8941	Sm with lg Y-valve	R669052	R200542	R2730	35mm
R8942	Sm with lg T-valve	R669052	R200542	R2700	35mm
R8944	Sm with med T-valve	R669052	R200543	R2600	35mm
R8946	Sm with sm T-valve	R669052	R200591	R1410A	22mm
<b>Pedicatric Mouth &amp; Nose Breathing Masks</b>					
R8951	Lg ped with med Y-valve	R669053	R200584	R2630	35mm
R8952	Lg ped with med T-valve	R669053	R200584	R2600	35mm
R8953	Lg ped with sm Y-valve	R669053	R200585	R1420A	22mm
R8954	Lg ped with sm T-valve	R669053	R200585	R1410A	22mm
R8963	Sm ped with sm Y-valve	R669072	R200585	R1420A	22mm
R8964	Sm ped with sm T-valve	R669072	R200585	R1410A	22mm
V870C	Clip for mask strap, 1 only				

Last update 3/2004